Rupture of Membranes Test

INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
GEBRAUCHSANWEISUNG
GEBRUIKSAANWIJZING
ISTRUZIONI PER L'USO
KULLANIM TALIMATLARI

ROM-1025

TABLE OF CONTENTS

EN	Instructions For Use	3
ES	Instrucciones de uso	9
DE	Gebrauchsanweisung	15
NL	Gebruiksaanwijzing	22
IT	Istruzioni per l'uso	28
TU	Kullanım Talimatları	35

INTENDED USE

The Clinical Innovations ROM Plus ES fetal membrane rupture test is a rapid, qualitative immunochromatographic test for the in vitro detection of amniotic fluid in vaginal secretions of pregnant women with signs and symptoms of ROM. The test detects AFP (alpha-fetoprotein) and IGFBP-1 (insulin-like growth factor-binding protein-1 also known as PP12 (placental protein 12)) from amniotic fluid in vaginal secretions. The test is for prescription use by health care professionals to aid in the detection of rupture of membranes (ROM) in pregnant women in conjunction with other signs and symptoms.

REAGENTS AND COMPONENTS

Each ROM Plus ES kit contains:

- One ROM Plus ES test strip
- One sterile polyester vaginal swab
- One vial with buffer solution (phosphate-buffered saline)
- Directions for use (on each pouch)
- The test strip contains the following antibodies: sheep antihuman IGFBP-1, goat antihuman AFP, mouse conjugated colloid gold antihuman IGFBP-1 & AFP.

TEST SUMMARY

ROM Plus ES is a self-contained test kit that provides qualitative results for rupture of fetal membranes and is obtained at point-of-care sites. A speculum is not needed to obtain ROM Plus ES results. The test is non-invasive, with only a simple vaginal swab sample required. The sample is collected by placing the swab into the vagina for 15 seconds. The swab is then mixed into a vial containing 600 μ L of buffer solution for 15 seconds and the ROM Plus ES test strip is placed into the buffer solution. The liquid moves chromatographically and unidirectionally towards the absorbent pad.

TEST PRINCIPLE

During migration, the sample reacts with mono/polyclonal antibodies bound to the test strip membrane. These antibodies are immunoreactive to a combination of proteins, IGFBP-1 and AFP, which are markers of amniotic fluid. As the membrane absorbs the liquid sample, a control line will appear, indicating proper function of the device. If the sample contains the IGFBP-1 and/or AFP proteins, it binds to the antibody of the test line, causing the test line to appear and indicating a positive result. If the sample does not contain the IGFBP-1 and/or AFP, only the control line will be visible, indicating a negative result.

STORAGE AND STABILITY

- Store the kit in a dry place at 4° to 37°C (40° to 99°F). Do not freeze.
- When stored in the foil pouch at the recommended temperature, the test is stable until the "Expiration" date on the pouch.

- Use ROM Plus ES within six (6) hours after opening foil pouch.
- Use ROM Plus ES within six (6) hours of collecting the vaginal swab sample and placing it into the buffer vial.

PRECAUTIONS, LIMITS AND WARNINGS

- ROM diagnoses should not be based on any single test.
- ROM Plus ES is for in vitro diagnostic use only.
- ROM Plus ES is for healthcare professional use only.
- Allow pouch containing ROM Plus ES to reach room temperature prior to utilizing the test.
- All instructions should be followed carefully for accurate results.
- Each ROM Plus ES test kit is single use and disposable and should not be reused.
- ROM Plus ES results are qualitative. No quantitative interpretations should be made.
- ROM Plus ES test kits will function properly with trace amounts of blood in the sample. Significant amounts of blood discharge may cause the test to malfunction and is not recommended.
- Safety precautions should be observed when collecting, handling, and disposing of test samples.
- Used test kits are biohazardous.
- Elevated fetal serum, urine, cord blood, and amniotic fluid as well as maternal serum levels of AFP have been reported in the literature in various developmental disorders such as neural-tube defects, hypothyroidism, autoimmune states, congenital heart defects, cystic fibrosis, etc. ROM Plus ES has not been evaluated for potential interference in these conditions.

Results should be used in conjunction with other clinical information available

METHOD OR PERFORMANCE

The ROM Plus ES assay was validated for the parameters of linearity, limit of detection, accuracy/reproducibility, sensitivity, specificity and cross-reactivity:

High Concentration (“High Dose Hook” effect) – for the ROM Plus upper-detection range, the IGFBP-1 and AFP were tested. Concentrations of IGFBP-1 were tested up to 400,000 ng/ml and AFP up to 200,000 ng/ml with a positive visual results. (Although lines may be lighter in presence of very high AFP and IGFBP-1 concentrations, any line is considered positive.)

The lowest limit of detection (LOD) is 5 ng/ml for IGFBP-1 and 150 ng/ml for AFP. (These LOD concentrations refer to those in samples of vaginal secretions before dilution with buffer.)

Reproducibility was tested on different days at six levels of amniotic fluid spiked

into a negative control. The assay was run on three lots of ROM Plus to determine the visual positive or negative results. Two low positives, two moderate positives and two high positives were run on three lots of ROM Plus on four different days. No difference in activity was observed.

To determine interference and cross-reactivity of the assay, Tylenol, aspirin, and five different bath products (KY Gel, Surgilube, Lever Soap, Noxema cream, Pert Shampoo), were spiked into the low positive control at a final concentration of 0.1% without visual loss of activity. The same bath products were spiked into the negative-matrix control and shown to be negative. In addition, human semen, urine and blood were spiked into the low positive at a 10% final concentration without loss of activity. Human semen, urine and blood were also spiked into the negative-control matrix and shown to be negative. IGFBP-1 assay does not cross-react with IGFBP-2, IGFBP-3, and IGFBP-4 based on Western Blot results. ROM Plus was shown to be negative when tested with specimens that were positive for bacterial vaginosis and sexually transmitted diseases. All samples were pH>4.5.

In a multi-center prospective, observational study, ROM Plus was shown to outperform the standard clinical assessment. Included in the study were healthy pregnant women, 18 years of age or older, between 15-42 weeks of gestation reporting signs or symptoms of rupture of membranes. Patients with known placenta previa and/or active vaginal bleeding were excluded from the study. Initial evaluation included both the standard clinical assessment for rupture of membranes (SSE: speculum exam for fluid pooling, ferning, and nitrazine test), as well the ROM Plus. Rupture of membranes was diagnosed if fluid was seen leaking from the cervical os, or suspected if two of the three conditions were present: pooling of fluid, positive nitrazine test, or ferning. The diagnosis of rupture of membranes was confirmed after review of the patient's clinical course following delivery.

Clinicians performing the SSE were blinded from the ROM Plus results and vice versa. The clinical standard of pooling/ferning/nitrazine has been shown in the literature to have sensitivities of 16-51% and specificities around 70%^{1,2,3} (when these tests are used individually) and thus are not ideal reference standards. As a result, the final diagnosis of ROM was based on subsequent review of the patient's chart and clinical course—a common methodology described in the literature for the detection of ruptured membranes.^{4,5}

Table 1: Clinical Study⁶ - ROM Plus vs. Clinical Assessment (pooling/ferning/nitrazine)

Of the 285 patients (15-42 weeks of gestation), the false positive rate for ROM Plus was 9% and the false negative rate was 0.5%, sensitivity 99%, specificity 91%, positive and negative predictive values of 95% and 99%. Conventional clinical evaluation's sensitivity was 85%, specificity 98%, positive and negative predictive values of 99% and 77%. Results in this study for ROM Plus in preterm were: sensitivity of 100%, specificity of 94%. Table 1 shows the ROM Plus test results (by gestational age) obtained in the clinical study. Six of the nine false positives were at term and only three of the nine were preterm. Data are

presented as sensitivity and specificity with confidence intervals for all gestational age group, term and preterm groups, less than 24 weeks' gestational age, between 24 and 34 weeks, and greater than 34 weeks gestational age.

GA (weeks)	N	Sensitivity (95% confidence interval)	Specificity (95% confidence interval)
ALL	285	99.5(97.0, 99.9)	90.7(88.3, 95.0)
>37	198	99.4(96.4, 99.9)	86.0(72.7, 93.4)
>37	87	100(89.6, 100)	94.4(84.9, 98.1)
>34	236	99.4(96.7, 99.9)	88.5(76.4, 93.8)
24-34	36	100(83.2, 100)	100(89.0, 100)
<24	13	100(43.9, 100)	80(49.0, 94.3)

QUALITY CONTROL

Each ROM Plus ES test contains Internal Quality Control, which verifies the integrity of the test procedure (operator), verifies proper assembly of the test strip (device integrity) and that the test has been appropriately stored (not exposed to extreme temperature) (environment).


The positive control, demonstrated by the appearance of the control line (C), assures that the reagent is functioning properly (device integrity), an adequate sample volume was applied to the cassette (operator) and adequate capillary migration (lateral flow) occurred (device integrity). The negative control is represented by the lack of color along the entire strip length assuring no inappropriate binding occurred (device integrity). If the control line appears and there is color along the strip the test is functioning properly and does not need to be repeated.

Because this test contains two levels (positive and negative) of internal quality control and is a qualitative, moderately complex test, it is the company's recommendation that external QC be completed for each new lot number or shipment of test material and when there's suspicion of compromised specimen or improper storage and as required by regulatory or accrediting agencies.

ROM Plus Fetal Membrane Rupture Test Quality Control kit for external controls are available from Clinical Innovations, LLC. and are intended to monitor the performance of Clinical Innovation's ROM Plus fetal membrane rupture test strip. The control test kits contain human IGFBP-1 and human AFP amniotic fluid proteins at concentrations of 20 ng/ml and 600 ng/ml. External negative controls that do not contain these proteins are also available in the Quality Control Kit. Always follow federal, state and local guidelines for quality control documentation.


TEST PROCEDURE

1 - Collect Sample


Remove ROM Plus ES contents from the packaging. Holding the buffer vial in an upright position, remove the shipping cap and set it aside. Remove the sterile swab from its package to collect a sample from the surface of the vagina. The tip of the swab should not touch anything prior to its insertion. Insert the swab tip into the vagina 2-3 inches (5-7 cm) deep. Withdraw the swab after a minimum of 15 seconds.

2 - Mix Sample and Run Test


Place the swab tip into the vial. Mix the swab in the buffer for at least 15 seconds, then remove and dispose of the swab.

3 - Add test strip to vial and set timer

Tear open the foil pouch and remove the ROM Plus test strip. Place the white end of the test strip (marked with arrows) into the vial with buffer.

4 - Read test result

Remove the test strip if two lines are clearly visible or after 10 minutes. Darkness of the stripes may vary. The test is still valid even if the lines are faint. Do not interpret test results based on darkness of the lines. If only a control line (top line) is visible, the test result is negative. If both the control line and test line (bottom line) are visible, the test result is positive. If no lines are visible the test result is invalid and should be repeated. A light visible line located in the test line region should be

considered positive; in addition, very high concentrations of proteins may result in a light test line. It is recommended to read the strip by 10 minutes after placing it into the vial.

EXPECTED RANGES

IGFBP-1 concentration in amniotic fluid as determined in literature studies is between 10,500 and 350,000 ng/ml and for AFP from 2,800 to 26,000 ng/ml and in serum from 55 to 242 U/ml (equivalent to 33 to 290 ng/ml). Concentrations of IGFBP-1 in amniotic fluid was observed to be 100 to 1000 times higher than that in maternal serum. These studies also showed that concentrations of these proteins in urine to be negligible.⁷⁻⁹

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Bibliography

1. Thomasino T, Levi C. Diagnosing PROM with Combination Monoclonal/Polyclonal Immunologic Protein Detection. ACOG 2011 Annual Clinical Meeting, Washington DC, Poster Session, May 2011.
2. Caughey AB, Robinson JN, and Norwitz, ER. Contemporary diagnosis and management of preterm premature rupture of membranes. *Rev Obstet Gynecol.* 2008;1(1):11-22.
3. Tagore S, Kwek K. Comparative analysis of insulin-like growth factor binding protein-1 (IGFBP-1), placental alpha-microglobulin-1 (PAMG-1) and nitrazine test to diagnose premature rupture of membranes in pregnancy. *J Perinat Med* 2010;38:609-12.
4. El-Messidi A, Cameron A. Diagnosis of Premature Rupture of Membranes: Inspiration From the Past and Insights for the Future. *J.Obstet Gynaecol Can* 2010;32(6):561-69.
5. Albayrak, M, Ozdemir A, Koc O, Ankarali H, Ozen O. "Comparison of the diagnostic efficacy of the two rapid bedside immunoassays and combined clinical conventional diagnosis in prelabour rupture of membranes," *Eur.J. of Obstet. & Gyn. and Reprod. Biology.* 158 (2011)179-182.
6. Cousins LM, Smok DP, Lovett SM, et al. AmniSure placental alpha microglobulin-1 rapid immunoassay versus standard diagnostic methods for detection of rupture of membranes. *Am J Perinatol.* 2005;22(6):317-320.
7. Lee SE, Park JS, Norwitz ER, et al. Measurement of placental alpha-microglobulin-1 in cervicovaginal discharge to diagnose rupture of membranes. *Obstet Gynecol.* 2007;109(3):634-640.
8. Lee SM, Lee J, Seong HS, Lee SE, et al. The clinical significance of a false positive AmniSure test in women with term labor with intact membranes. *J. Matern Fetal Neonatal Med* 2009 Apr; 22(4): 305-310.
9. Chen FC, Dudenhausen JW. Comparison of two rapid strip tests based on IGFBP-1 and PAMG-1 for the detection of amniotic fluid. *Am J Perinatol.* 2008;25(4):243-246.

Clinical Innovations
747 West 4170 South
Murray, Utah USA 84123
P: 801-268-8200
Toll Free: 888-268-6222
clinicalinnovations.com

Clinical Innovations Europe, Ltd.
6-9 The Square | Stockley Park,
Heathrow UB11 1FW | United Kingdom
Tel +44 (0) 333 370 4408

LATEX FREE, DEHP FREE | Made in the USA.
Patents Pending | Copyright 05-2016, Clinical
Innovations, LLC. | All rights reserved.

**CLINICAL
INNOVATIONS**

for MOM. for BABY. for LIFE.

Copyright 2016, Clinical Innovations, LLC.
All rights reserved. P/N 050-0742 REV C