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## Cleaning and Disinfection of Laborie UniTip Catheters

The instructions for use are authoritative for a correct cleaning / disinfection. This document is not part of the device's labelling. The latest version of this document is available on our website. Please check it with the appropriate frequency to be kept informed on the latest updates.

In addition to the recommended cleaning and disinfection solutions in the UniTip *Instructions for Use*, UniTip Catheters may be cleaned and disinfected with the listed solutions/systems **without damaging the UniTip Catheters (material compatibility and performance)**. Laborie makes no claims on the efficacy of cleaning and disinfection solutions listed in this document, as this is done by the manufacturer of each solution.

The recommended immersion time of the authoritative instruction for use must be maintained.

It's imperative to thoroughly rinse the UniTip Catheters after each cleaning and disinfection, regardless of the instructions of the solution manufacturer. This also applies to all Wipes Systems.

Because some automated reprocessors cannot generate the pressure required to force fluid through the lumen, manual cleaning and disinfection is the safest method to effectively reprocess this area of a lumened device.

After manual cleaning and disinfection of the lumen, mechanical reprocessing of the catheter is subsequently possible.

**Exception:** Mechanical reprocessing of lumened catheters using the Tristel Stella B-LP (disinfection cycle only) and Stella C-LP (cleaning and disinfection cycle) does not require manual cleaning and disinfection.

The USER is responsible for cleaning / disinfection / sterilization of the UniTip catheter. Take care that only proven procedures and solutions are used for cleaning / disinfection / sterilization, and that the parameters always comply with the validation.

Additionally, please pay attention to the laws / regulations / guidelines of your country as well as the hygienic instructions of your hospital.

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(Specifically tested and recommended for UniTip High-Resolution Catheters)

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### Mechanical Procedure (Low Temperature-Cycle <60°C)

#### UniTip High-Resolution Catheter with connector IPX5 / IPX7

(Unless otherwise noted, lumen must be reprocessed manually.)

**Hersteller / Manufacturer: Wassenburg Medical Devices B.V.; Gerät / Device: WD 440/WD 440 PT**

Reinigungslösung / Cleaning Solution: neodisher <sup>®</sup> MediClean forte	Dr. Weigert
Desinfektionslösung / Disinfection Solution: neodisher <sup>®</sup> Septo DN	Dr. Weigert

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**Hersteller / Manufacturer: Tristel ; Gerät/Device: Stella C-LP\* (Cleaning and Disinfection cycle)**

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\*Does not require manual lumen cleaning and disinfection