



Gynecology/Urology Procedure Chair

SONESTA **6202**



The 6202 procedure chair... the ultimate in positioning and comfort for you and your patient.



STILLE

The word "STILLE" is written in a large, bold, grey, sans-serif font. Above the text is a large, thick, red curved graphic element that arches over the letters.



Gynecology/Urology Procedure Chair

SONESTA 6202



The Stille 6202... Great Positioning For:

- Urodynamics
- General Examinations
- Cystoscopy
- Ultrasound
- Biofeedback
- Vasectomy
- OB/GYN Procedures

Engineered To Be Flexible:

- Locking mechanism for precise positioning
- Three motors for height, seat and back adjustments
- Easy positioning for elderly or physically challenged patients
- Wide range of adjustment, from a low of 16" (40cm) to 39.4" (100cm) high
- Engineered ergonomics. Infinitely adjustable (see specifications) for physician and patient comfort.
- Lightweight and mobile easy to move to other rooms; unmatched access for floor cleaning and disinfecting.
- Tough and reliable, able to support up to 375 lbs - 170kg for years of service in the most demanding practices.

Standard Accessories:

- 3 motor foot control 524-6357
- 1 seat cushion with 170mm cutout 524-6366-3
- Paper roll holder 524-6326
- 2 fold up footrests 524-6331
- 2 heavy duty Rissler leg rests 524-6336
- Drainage Receptacle 524-6325-1
- 2 hand grips/arm rests 524-6339-2
- 2 leg rest holders 524-6338
- 1 neck pillow 524-6346
- 1 full seat (w/o cut out) 524-6365

The 6202 chair is built and engineered to be the ultimate in positioning for your practice.

The 6202 offers the ultimate in positioning flexibility.

Packaged with our most popular accessories, the 6202 is fully equipped for a myriad of exam and treatment procedures. Three powerful (yet quiet) motors allow fast & precise positioning with minimal effort.

SPECIFICATIONS:

| | | |
|-------------------------------|--------|---------|
| LENGTH WITH EXTENSION CUSHION | 190cm | 75" |
| LENGTH WITHOUT | 143cm | 56" |
| WIDTH | 62cm | 25" |
| HEIGHT (MIN) | 40cm | 16" |
| HEIGHT (MAX) | 95cm | 37" |
| LIFT CAPACITY | 170kg | 375lbs |
| WHEEL DIAMETER (LOCKING) | 10cm | 4" |
| NET WEIGHT | <125kg | <276lbs |

ELECTRICAL:

120/230v 24v DC Output
 Meets EMC (Electromagnetic Compatibility) for IEC 60601

All Products carry the CE mark.

DISTRIBUTED IN U.S. BY:

Laborie Medical

400 Avenue D, Suite 10
 Williston, VT 05495
 Phone: +1 800.522.6743
 Fax: +1 802.878.1122
 Web: www.laborie.com

CONTACT IN CANADA:

Laborie Medical

2101 Boulevard Lapinière
 Brossard, Québec J4W 1L7
 Phone: +1 888.522.6743
 Fax: +1 450.671.7182
 Web: www.laborie.com

