

delphis KT →



DELPHIS FAMILY // QUALITY - INNOVATION - SERVICE - SATISFACTION

LABORIE
urodynamics

TECHNICAL SPECIFICATIONS



SOFTWARE SPECIFICATIONS & OPTIONS

PATIENT MANAGEMENT SOFTWARE

Patient management software enables a consolidation of UDS™ data and ICS approved questionnaires, database queries and customized report generation as well as other time-saving features.

INTUITIVE WINDOWS BASED SOFTWARE

Complete and versatile Windows based software makes your everyday urodynamics testing quick, reliable and simple.

VBN MODELIZED ANALYSIS

VBN exclusive modeling software provides visual comparison of flow curves.

EXTENSIVE TESTING CAPABILITIES

All ICS standardized tests can be performed with confidence: **Uroflowmetry, Cystometry (CMG), Urethral Pressure Profile (UPP), Leak Point Pressure (LPP)** and complete micturition studies with **EMG**.

INFORMATION TECHNOLOGY

I.T solutions enable integration of UDS findings with information systems in order to facilitate the patient care management process.

PRODUCT SPECIFICATIONS & OPTIONS

RECORDING CHANNELS

- Two (2) pressure transducers
- Built-in mini pump
- Uroflow transducer with commode chair
- EMG unit
- UPP puller**
- Third pressure transducer**
- Fourth pressure transducer**
- Anorectal manometry**

PROCESSOR / PRINTER

Laptop with color deskjet printer.

MODULAR CART

A customizable modular unit contained in an attractive, ergonomically-sound cart. An ideal work station featuring in-office portability.

** Optional

UDS™ PERIPHERALS & ACCESSORIES

UROFLOW TRANSDUCER

With its magnetic dish, the Uroflowmeter is sensitive enough to detect low flow rates as well as tolerate high stress.

PRESSURE TRANSDUCERS

Compatible with fluid-filled, air-charged and electronic catheters to monitor vesical, abdominal and urethral pressures. Wide range of catheters including the Laborie™, Unisensor™ & T-DOC™ catheters.

MINI PUMP

Integrated pump provides smooth operation with high accuracy and variable filling rates. Volume and pressure limit warnings for added safety & support.

UPP

Software controlled puller for greater precision. Automatic UPP event annotations, UPP summary including peak PURA & peak PCLO calculations with up to 8 pull capacity.**

ELECTROMYOGRAPHY

High-fidelity integrated EMG is compatible with several styles of surface and needle electrodes. Flexible software allows for customizable EMG displays.

** Optional

LABORIE

urodynamics



USA

LABORIE MEDICAL TECHNOLOGIES CORP.
400 Avenue D, Suite 10
Williston, Vermont
05495-7828
Tel.: +1 802 857 1300
800 522 6743

CANADA

LABORIE MEDICAL TECHNOLOGIES INC.
2101 Boulevard Lapinière
Brossard, Québec
J4W 1L7
Tel.: +1 450 671 5901
888 522 6743

INTERNATIONAL

LABORIE MEDICAL TECHNOLOGIES INC.
6415 Northwest Drive, Unit 10
Mississauga, Ontario
Canada L4V 1X1
Tel.: +1 802 857 1314

EUROPE

LABORIE MEDICAL TECH. N.V.
Technologiepark, Zwijnaarde 3
9052, Gent, Belgium
Tel.: +32 9 222 39 85

CHINA

LABORIE CHINA
Rm.31B BLDG No.3
Century Classic, No. 45
Xiao Guan Bei Li, Beijing
100029
Tel.: +861 084 898 891

LABORIE REPRESENTATIVE

LABORIE

For product demonstration,
please contact us at:
marketing@laborie.com
and visit:
www.laborie.com

ISO 13485:2003