



# UROCAP™ IV

INNOVATING THE WAY URODYNAMICS IS DELIVERED



Innovation for Health

# UROCAP™ IV

## UROFLOWMETRY REFINED

UROCAP IV is a reliable, standalone, cost effective and easy to use Uroflowmetry solutions from LABORIE.

As a wireless, battery-powered, water-protected device, it is capable of providing patient privacy and functional versatility; simply place it where you want and it will send data wirelessly to your computer.

The UROCAP IV software, combined with iList Reporter software, is easy to use and allows you the ability to review and generate a variety of scripted or data reports with a few clicks of a mouse! With one extra click, the reports can be printed, sent to MS Word, PDF writer or automatically emailed to referring physicians. Our latest technology allows you to simplify your clinical workflow and increase patient throughput without strain on your team or compromise on quality of care.

As your solution provider for quality Uroflowmetry, LABORIE's UROCAP IV is designed with your work flow in mind.



## TECHNICAL SPECIFICATIONS

Size & Weight	L = 6.3" (16.0cm)    W= 6.3" (16.0cm) H = 2.6" (6.6cm)    Weight = 0.7lbs (0.3kg)	
Measurement Channels	Volume: Range = 0 to 1200ml Accuracy = +/- 2% Sampling Rate = 20Hz	Flow: Range = 0 to 50mL/s Accuracy = +/- 5% Sampling Rate = 20Hz
Indicator Lights	Bluetooth, Status & Battery	
Charger	5 VDC, 3.0A	


## OPTIONAL ACCESSORIES



Commode Chair



Uroflow Stand



**USA:**  
Tel.: 800 522 6743  
Email: [usmarketing@laborie.com](mailto:usmarketing@laborie.com)

**EUROPE/INTERNATIONAL:**  
Tel.: +31 53 4803700  
Email: [info@laborie.com](mailto:info@laborie.com)

**HEADQUARTERS:**  
Tel.: +1 905 612 1170  
Email: [marketing@laborie.com](mailto:marketing@laborie.com)



[www.laborie.com](http://www.laborie.com)