



PORTASCAN™ 3D

PORTABLE 3D BLADDER SCANNER



Innovation for Health

AN IMAGE YOU CAN RELY ON

INTRODUCING PORTASCAN™ 3D

FEATURES AND BENEFITS

- Automatic 3D Bladder Scanner for fast and accurate bladder volumes
- Easy-to-use targeting and real-time scanning modes, ideal for non-specialists or experienced users
- Calibration not required, reducing time and cost of annual maintenance
- Touchscreen interface enhances procedure simplicity
- Light weight ultra-portable design
- Easy point-and-freeze operation provides results within 5 seconds
- High capacity internal memory for saving patient data
- Simple-to-use, ergonomically designed probe

SPECIFICATIONS

Procedure Type	Bladder volume and post void residual scanning (PVR)
Monitor Type	LCD Color Touch Screen
Monitor Size	6 in
Printer	Thermal
Battery Capacity	Up to 4 hours
PC Connectivity	USB
Dimensions	14 cm x 32.5 cm x 15.5 cm (5.5 in x 13 x 6 in)
Weight	2 kg (4 lbs.)
Storage	100 scans



USA:
Tel.: 800 522 6743
Email: usmarketing@laborie.com

EUROPE/INTERNATIONAL:
Tel.: +31 53 4803700
Email: info@laborie.com

HEADQUARTERS:
Tel.: +1 905 612 1170
Email: marketing@laborie.com



www.laborie.com