



DUAL-PURPOSE EMG AND STIMULATION PROBES

PELVIC FLOOR REHABILITATION AND INCONTINENCE TREATMENT



Innovation for Health





COMPREHENSIVE RANGE OF PROBES FOR PELVIC FLOOR REHABILITATION


PROBES YOU HAVE TRUSTED FOR MORE THAN 20 YEARS

LABORIE is a global market leader in pelvic health diagnostic and treatment, offering vaginal and anal rehabilitation/biofeedback probes. These probes provide comfortable muscle stimulation and accurate biofeedback for male and female patients. LABORIE's large selection of probe shapes and sizes, allows for use with most commonly available equipment.

SELECT THE PERFECT CONFIGURATION FOR EACH PATIENT

Choose your desired probe type to find the correct part number.

PROBE DESCRIPTIONS	PROBE SHAPES	
VAGINAL CLASSIC PROBE <i>For common use</i> Length is 110mm Diameter is 25mm		P/N: PRB450
VAGINAL TAMPON PROBE <i>For short vaginas</i> Length is 50mm Diameter is 20mm		P/N: PRB425
VAGINAL ATROPHIC PROBE <i>For narrow vaginas</i> Length is 110mm Diameter is 19mm		P/N: PRB085
ANAL SHORT PROBE <i>For common use</i> Length is 83mm Diameter is 14mm		P/N: PRB055



USA:
 Tel.: +1 802 857 1300
 800 522 6743
 Email: usmarketing@laborie.com

EUROPE/INTERNATIONAL:
 Tel.: +31 53 4803700
 Email: info@laborie.com

HEADQUARTERS:
 Tel.: +1 905 612 1170
 Email: marketing@laborie.com



ANNIVERSARY

www.laborie.com