



PrimeSight™ UNITY Cart

ENDOSCOPY CART



Innovation for Health

PRIMESIGHT™ UNITY ENDOSCOPY CART

KEY FEATURES

The PrimeSight™ UNITY Cart is designed to store the necessary equipment and instruments required for performing endoscopic procedures using the PrimeSight Cystoscopy system. The cart was designed with a small footprint to facilitate easy movement from room to room and includes the following design features:


- Tabletop height adjusts with the press of a lever
- Embedded antimicrobial powder coat finish
- Integrated design features include scope storage, EndoSheath loading platform and EndoSheath storage baskets
- Lockable caster wheels for stationary security
- 4-outlet medical grade power strip
- Easy grip handles
- Height adjustable IV/Saline/Water Pole
- Optional accessories include printer, IV pole, cart cover, keyboard, endoscope storage shield, and foot pedal.

DIMENSIONS (INCLUDING WHEEL BASE)

Width	671mm
Height (without IV pole)	942mm; adjusts to 1197mm
Height (with IV pole)	<1696mm; adjusts to 2256mm
Depth	611mm
Caster/wheel size	>125mm
Shelf size	428mm (W) x 328mm (D)
Drawer size	352mm (W) x 235mm (D) x 139mm (H)
Wire basket size	140mm (W) x 94mm (D)
Weight	<38.5kg (85 lbs)

CATALOG NUMBER	DESCRIPTION
07-9022	PrimeSight UNITY cart
CART ACCESSORIES	
07-9012	Keyboard
07-9013	Foot pedal
07-9024	IV pole, adjustable (795mm extends to 1100mm)
07-9025	Endoscope storage shield
07-9018	HP laser printer
07-9020	Video printer Sony UP-D25MD
07-9021	Video printer Sony UP-DR80MB
07-9026	Cart cover
07-9027	PrimeSight UNITY strap

PrimeSight is a trademark and EndoSheath is a registered trademark. UNITY Cart is manufactured by Cogentix Medical. ©2019 LABORIE. All rights reserved.



USA:
Tel.: +1 802 857 1300
800 522 6743
Email: usmarketing@laborie.com



www.laborie.com