

EndoTNeedle™

Endoscopic Injection Needle

GI Supply®
Specialty Endoscopic Products



TorturLESS Path

A tough needle, at a fair price.

EndoTNeedle was engineered to balance performance and value while delivering a highly reliable needle for your most tortuous cases. Designed with quality materials, features from premium needles, and a handle your staff will love. The EndoTNeedle is your new everyday needle at a price you are willing to pay.

Click It

- Your nurses will love the ergonomic handle for single-hand use
- Thumb activated needle extension and spring loaded retraction

Stick It

- Reliable performance in retroflex and tortuous paths
- Stainless steel sheath prevents kinking
- Smooth round hub which aids in the closure and blockage needed to stop GI bleed cases
- Ultra sharp needle



Order Number	Minimum Working Channel	Working Length	Needle Extension	Needle Gauge	Units/Box
ET1300	2.8mm	160cm	4mm	25	10
ET1301	2.8mm	230cm	4mm	25	10
ET1302	2.8mm	160cm	5mm	25	10
ET1303	2.8mm	230cm	5mm	25	10
ET1320	2.8mm	160cm	4mm	23	10
ET1321	2.8mm	230cm	4mm	23	10
ET1322	2.8mm	160cm	5mm	23	10
ET1323	2.8mm	230cm	5mm	23	10

The perfect fit, in a kit!

Spot® Ex and EndoTNeedle go together perfectly. Now you can order both with one item number. The kit includes a box of Spot Ex and a box of the EndoTNeedle.

Item number GIS-4501. Add EndoTNeedle to your next GI Supply order.