

FLUID-FILLED CATHETERS

ACCURATE – PRACTICAL – COMPREHENSIVE



Innovation for Health

SINGLE-USE URODYNAMIC CATHETERS


PRECISE DIAGNOSTIC CONFIDENCE

STRONG AND FLEXIBLE. Assembled with thermoplastic material that softens over time for easier insertion and patient comfort. Unique wall material eliminates crimping and folding of the catheter during insertion.






UNIQUE DESIGN. The molding construction has a precise design to allow the correct flow without turbulence between the catheter body and the infusion line to prevent pressure artifacts.


USER FRIENDLY. Centimeter graduations to ease positioning of catheter. Luer lock connectors are compatible with most urodynamic systems and color coded for fast procedure setup.

COMPLETE RANGE. Extensive family of options designed for different patients and procedures. Choose the French size, tip style and number of lumens to suit your needs.



*Please refer to the catalog for complete range of French sizes, tips and options

| CATHETER PORT COLOR | DESCRIPTION | CONNECTION |
|--|----------------|------------------|
| Clear  | Infusion line | Pump Tubing |
| Blue  | Vesical line | Pves |
| Red  | Urethral line | Pura |
| White  | Abdominal line | Pabd |
| Yellow  | Perfusion line | Perfusion Tubing |



USA:
Tel.: 800 522 6743
Email: usmarketing@laborie.com

EUROPE/INTERNATIONAL:
Tel.: +31 53 4803700
Email: info@laborie.com

HEADQUARTERS:
Tel.: +1 905 612 1170
Email: marketing@laborie.com



www.laborie.com