



LABORIE PRODUCTS HELP GUIDE THE PATHWAY TO BETTER UROLOGICAL HEALTH

LABORIE helps providers along the care pathway with a comprehensive suite of products and services and offers options for the diagnosis, treatment, and ongoing support of their patients.



The NXT™ Pro is part of the Laborie Intelligent Cath-ID System

Supporting Patient Safety and Enhancing Urodynamics Quality

- Supports best practices in single patient use of consumables to minimize risk of cross-infection and/or cross-transmission
- Catheter compatibility for quality urodynamic studies
- Supports traceability, removes need for manual entry, reducing errors while automating workflow


WE ARE HERE TO SUPPORT YOU

Visit us at: <https://www.laborie.com/product/NXT-range>

Talk to your Laborie Territory Manager about how NXT™ can help your practice.

Some Laborie products may not be available in certain regions. Please contact your local Laborie Sales Representative to learn more about the products available to you.
© 2023 Laborie. All rights reserved.

Laborie ¹	
USA	EUROPE/INTERNATIONAL
T +1 800 522 6743	T +31 53 4803700
E marketing@laborie.com	E marketing@laborie.com



Learn more

1 01394 - 2

Laborie¹
FOR DIGNITY. FOR LIFE.



NXT™ Pro

Defining the New Standard in Urodynamics

UROLOGY & UROGYNECOLOGY
NXT™ is the new standard for modern Urodynamics to improve physician confidence and patient experience.



NXT™ Systems
The New Standard in Urodynamics

NXT™ Pro System

NXT™ Pro is an advanced Urodynamics System offering a complete solution, including full video, HL7 connectivity and an expandable platform for future software functionality.

Infusion Hook:
can be rotated 180 degrees from front to back

Flexible Hoop:
Choose right or left-handed configurations

Roams:
NXT Pro comes with two roams to allow interchange of air and fluid set-ups efficiently

Consumable Scanning:
for documentation and traceability

Water-resistant Components:
for easy cleaning and wipedown in between studies

Catheters:
Intuitive design for easy setup and procedure start.

Touchscreen:
Navigation without a keyboard and mouse

Remote Control:
Mark key events from across the room

Storage:
A drawer is integrated into the top shelf for organization of additional items

Adjustable Spine:
Move to your ideal height

Small Footprint:
Easy to roll between rooms

Synergy Software

URODYNAMICS, YOUR WAY.

Customizable Workflows and Reports:



Configurable workflow and reporting tools address the unique needs of academic facilities, hospitals and clinics.

Enhanced User Experience:

Intuitive design, workflow guidance and seamless interface puts NXT™ at the forefront of modern urodynamics.

Next Generation Decision Support and Error Prevention:

Guided workflow helps in training with new staff and ensuring test consistency.



The **NXT Urocap** comes wired or wireless and designed for seamless integration within the NXT platform.



The **NXT Pump module** features state of the art technology, a wide range of flow rates and multiple calibration line storage.



The **NXT all-in-one computer** offers a novel touchscreen experience.



A **wireless capable printer** included for convenient document creation and immediate review with patients.



The **NXT Roam and PIM** is a sophisticated module that delivers pressure and EMG readings in a user-friendly format.



The **UPP Puller** boasts Laborie's smallest footprint and completely wireless battery operation. Available on a stand-alone rolling floor stand.