

# EMPOWER YOUR DATA

## SOFTWARE & DATA MANAGEMENT

### CUSTOMIZABLE PROCEDURES AND REPORTS

- Comprehensive range of preset options can be modified as desired

### TEST PLAYBACK IN REAL-TIME OR ACCELERATED MODES

- Easy review and teaching with LABORIE systems

### HOT KEYS & VOICE ANNOTATIONS TO GUIDE USER THROUGH STUDY

- Helps clinician follow standardized workflow

### TRUE EMR & HL7 INTEGRATION

- Integrate your data seamlessly through the proven i-LIST™ interface
- Easily manage and integrate data

### GOLD STANDARD REPORTING


- Choose from an array of comprehensive reports with just one click to minimize time spent on analysis



The Goby is part of the Laborie Intelligent Cath-ID System

### Supporting Patient Safety and Enhancing Urodynamics Quality

- Supports best practices in single patient use of consumables to minimize risk of cross-infection and/or cross-transmission
- Catheter compatibility for quality urodynamic studies
- Supports traceability, removes need for manual entry, reducing errors while automating workflow



Learn more

Goby-ARUEN-US-101875-1

**Laborie**<sup>7</sup>  
FOR DIGNITY. FOR LIFE.



## Goby

Compact · Portable · Wireless

### UROLOGY & UROGYNECOLOGY

The Goby is portable, modular and wireless device offers a unique combination of features to enhance your practice.



# SMALL IS BEAUTIFUL

## FEATURES AND BENEFITS

As the world’s most compact urodynamics system, Goby is LABORIE’s market-leading machine designed specifically for clinic and mobile use. This portable, modular and wireless device offers a unique combination of features to enhance your practice.

### FULL URODYNAMICS TEST CAPABILITY

- Perform all standard clinic procedures including Uroflow, Cystometrogram (CMG), Pressure Flow Study, Urethral Pressure Profiles (UPP), Stress Tests & Leak Point Pressures

### WIRELESS AND MODULAR DESIGN

- Maximize convenience for both clinicians and patients with the combination of battery power and Bluetooth® technology

### VARIABLE MOUNTING OPTIONS

- Choose the Goby set-up to suit your clinical needs with Cart, Procedure Chair, IV Pole, Tabletop or Travel Case options available

### OPTIONAL ACCESSORIES

- Enhance the system to make your workflow easier with a UPP Puller, Infusion Transducer, Commode Chair, Uroflow Stand or Remote Control

### ADDITIONAL APPLICATIONS

- Specialized options for different procedures including Clinic Urodynamics, Mobile Use, Goby EQ (EMG & Uroflow) and ARM (Anorectal Manometry)



UDS ROAM

GOBY HUB

UROCAP IV

CONFIGURE FOR YOUR CLINICAL NEEDS

## COMPONENTS

### UDS Roam

- Wireless Catheter Interface
- Records up to 4-Pressure & 1 EMG.
- Belt & Lanyard mountable to keep cables out of the way.

### Goby Hub

- Wireless Urodynamics Processor
- Precision Infusion Pump and Roam Charging Station.

### Urocap IV

- Wireless Uroflowmeter
- Remote connection offers patient privacy.
- Fourth generation water-resistant device with high quality recording.



Goby CT  
Cart



Goby IV  
IV pole



Goby SC  
Sonesta Chair



Goby TP  
Table Top



Goby EQ  
Flow & EMG option



ARM  
4-pressure Anorectal  
Manometry option



Flowstand  
Support Urocap  
IV flowmeter for  
standing use



Portable  
Backpack option