

# Eclipse™ System






for bowel control in women

UROLOGY & UROGYNECOLOGY







The Eclipse System is an innovative nonsurgical therapy that offers immediate results for women experiencing loss of bowel control.

## Ordering Information



|   |   |
|---|---|
|    | <p><b>Sizer</b></p> <p>Reusable sizing tool, representing the Insert base. Intended to aid in selecting appropriate Trial Insert sizes. Sizers are only to be used in the clinical setting and are not intended to be taken home.</p> <p>SZ02-48RD, SZ02-51RD, SZ02-57RD, SZ02-64RD, SZ02-70RD, SZ02-76EP</p>   |
|   | <p><b>Trial Insert</b></p> <p>Single-patient-use vaginal Insert intended for in-office use during fitting and evaluation process. The Trial Inserts are identifiable by their white color. (One Pump included per system and patient - patient should use same Pump with both Trial and Eclipse Inserts.)</p> <p><b>Small Balloon:</b> TR02-48RDS, TR02-51RDS, TR02-57RDS</p> <p><b>Regular Balloon:</b> TR02-48RDR, TR02-51RDR, TR02-57RDR, TR02-64RDR, TR02-70RDR, TR02-76EPR</p>                                     |
|  | <p><b>Eclipse Insert</b></p> <p>Single-patient-use vaginal Insert Intended for longer-term use (up to one year). The Eclipse Inserts are identified by their indigo color. (One Pump included per system and patient - patient should use same Pump with both Trial and Eclipse Inserts.) A new Eclipse device should be prescribed annually.</p> <p><b>Small Balloon:</b> EC02-48RDS, EC02-51RDS, EC02-57RDS</p> <p><b>Regular Balloon:</b> EC02-48RDR, EC02-51RDR, EC02-57RDR, EC02-64RDR, EC02-70RDR, EC02-76EPR</p> |
|  | <p><b>Pump</b></p> <p>The Pump is used by the patient for inflating and deflating the Insert Balloons. The Pump is packaged with the medium Regulator and a Patient Instructions For Use Booklet. (One Pump included with system and provided separately.)</p> <p>PP02</p>  |
|  | <p><b>Valve</b> (5 each per pack)</p> <p>Replacement Eclipse Valve, in the case of damage or irreversible soiling of the Valve. (Each Eclipse device includes one pre-installed Valve.)</p> <p>VV04</p>   |

## Ordering Information (continued)

|   |   |   |   |
|---|---|---|---|
|                                  | <p><b>Regulator</b> (5 each per pack)</p> <p>The medium Regulator is provided with the Pump. The low or high Regulator can be used to provide lower or higher maximum balloon pressure levels.</p> <p>RR02-LOW, RR02-MED, RR02-HIGH</p>   |   |   |
|                                  | <p><b>Extension Tube</b> (5 each per pack)</p> <p>Additional tubing to be used as needed to customize the Tube length insert.</p> <p>TT02</p>   |   |   |
|                                  | <p><b>Utility Cart</b></p> <p>4-drawer storage Cart to hold Eclipse products and patient materials.</p> <p>ECLCRT01</p>   |   |   |
|                                  | <p><b>Quick Start Kit (QSK) Options (one time purchase)</b></p>   |   |   |
|   | <p><b>Base Kit Only (No Eclipse Inserts)</b></p> <p>Includes:<br/>                     Sizers (11 Sizers – 2 each except 1 of the largest)<br/>                     Trial Inserts (14 Inserts – all 9 sizes plus extra of common sizes)<br/>                     Pump (3 Pumps)<br/>                     Accessories (5 each of the following:<br/>                     Valves, Extension Tubes, Regulator – low pressure, Regulator – high pressure)<br/>                     Utility Cart</p> | <p><b>5-Pack QSK</b></p> <p>Includes:<br/> <b>BASE KIT plus</b><br/> <b>5 Eclipse Inserts</b> (one each of the most common sizes)<br/>                     • 1 each: EC02-48RDR, EC02-51RDR, EC02-57RDR, EC02-64RDR, EC02-70RDR</p> | <p><b>14-Pack QSK*</b></p> <p>Includes:<br/> <b>BASE KIT plus</b><br/> <b>14 Eclipse Inserts</b> (all 9 sizes plus extra of the common sizes)<br/>                     • 1 each: EC02-48RDS, EC02-51RDS, EC02-57RDS, EC02-76EPR<br/>                     • 2 each: EC02-48RDR, EC02-51RDR, EC02-57RDR, EC02-64RDR, EC02-70RDR</p> |
|   | <p>ECLBDL01</p>   | <p>ECLBDL04</p>   | <p>ECLBDL03</p>   |
| <p><i>*Consignment option available for 14-Pack QSK - no upfront cost, but higher per unit purchase price</i></p> |   |   |   |
|                                | <p><b>5-Pack of Eclipse Inserts and Pumps</b> (special re-order bundle discount)</p> <p>Re-order Eclipse Inserts in quantities of 5 to receive bundle discount<br/>                     Account selects sizes of Eclipse Inserts for 5-pack re-order bundle (includes 5 Pumps (PP02))</p> <p>ECLBDL05</p>   |   |   |
|                                | <p><b>Eclipse Evaluation Program</b></p> <p>Laborie representative delivered product and case support to evaluate Eclipse System on a limited number of patients. Pricing pre-defined for the Trial Insert fitting and Eclipse Insert fitting inclusive of necessary accessories.</p>   |   |   |

To incorporate the Eclipse™ System into your practice, please visit us at [www.laborie.com/product/eclipse](http://www.laborie.com/product/eclipse) or scan QR code below.

Please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

©2023 Laborie. All rights reserved.

**Laborie**

USA

T +1 802 857 1300  
 +1 800 522 6743

E [usmarketing@laborie.com](mailto:usmarketing@laborie.com)  
 W [laborie.com](http://laborie.com)



Learn more