



LABORIE PRODUCTS HELP GUIDE THE PATHWAY TO BETTER UROLOGICAL HEALTH

LABORIE helps providers along the care pathway with a comprehensive suite of products and services and offers options for the diagnosis, treatment, and ongoing support of their patients.



The NXT™ Evo is part of the Laborie Intelligent Cath-ID System

Supporting Patient Safety and Enhancing Urodynamics Quality

- Supports best practices in single patient use of consumables to minimize risk of cross-infection and/or cross-transmission
- Catheter compatibility for quality urodynamic studies
- Supports traceability, removes need for manual entry, reducing errors while automating workflow

WE ARE HERE TO SUPPORT YOU

Visit us at: <https://www.laborie.com/product/NXT-range>

Talk to your Laborie Territory Manager about how NXT™ can help your practice.

Some Laborie products may not be available in certain regions. Please contact your local Laborie Sales Representative to learn more about the products available to you.
© 2023 Laborie. All rights reserved.


Laborie

USA

T +1 800 522 6743
E marketing@laborie.com

EUROPE/INTERNATIONAL

T +31 53 4803700
E marketing@laborie.com



Learn more

108177-1

Laborie
FOR DIGNITY. FOR LIFE.



NXT™ Evo

Defining the New Standard in Urodynamics

UROLOGY & UROGYNECOLOGY
NXT™ is the new standard for modern Urodynamics to improve physician confidence and patient experience.



NXT™ Systems

The New Standard in Urodynamics

NXT™ Evo System

NXT™ Evo includes many advanced urodynamics capabilities, while offering efficiency through user-centric design.

Remote Control:
Mark key events across the room.

Docking Station:
Intuitive single roam docking mount for storage & charging. Able to mount on left or right side.

Aesthetic Design:
Clean aesthetic design with electrical protection via LIT (Line Isolation Transformer).



Infusion Hook:
Left or right mounted infusion hook to conduct measurements during studies.


All-In-One Computer with Touchscreen:
Enables Navigation without a keyboard and mouse.

Portable Pump Module:
Enables consumable scanning and traceability.

Adjustable Shelf:
Manually adjustable shelf unlocked via turn knob.

Water-resistant components:
Allow for easy cleaning and wipe down between studies.

Catheters:
Intuitive catheter design for easy setup and procedure start.



Synergy Software

URODYNAMICS, YOUR WAY.

Customizable Workflows and Reports:



Configurable workflow and reporting tools address the unique needs of academic facilities, hospitals and clinics.

Enhanced User Experience:

Intuitive design, workflow guidance and seamless interface puts NXT™ at the forefront of modern urodynamics.

Next Generation Decision Support and Error Prevention:

Guided workflow helps in training with new staff and ensuring test consistency.



The **NXT Urocap** comes wired or wireless and designed for seamless integration within the NXT platform.

The **NXT Pump module** features state of the art technology, a wide range of flow rates and multiple calibration line storage.

The **NXT all-in-one computer** offers a novel touchscreen experience.

A **wireless capable printer** included for convenient document creation and immediate review with patients.

The **NXT Roam** and **PIM** is a sophisticated module that delivers pressure and EMG readings in a user-friendly format.

The **UPP Puller** boasts Laborie's smallest footprint and completely wireless battery operation. Available on a stand-alone rolling floor stand.