



LifeBubble[™]

Umbilical Catheter Securement System

NEONATAL

LifeBubble brings a new standard in the securement and protection of umbilical cord catheters.

Protecting the
lifelines of our most
vulnerable patients



Small




Large



LifeBubble[™]
Umbilical Catheter
Securement System

A Simple Way to Secure

LifeBubble is designed to secure umbilical catheters and protect the catheter insertion site for neonates in intensive care.



- **Reduces Risk of Catheter Migration and Early Discontinuation** with stronger securement^{1,2}
 - **Protects the Insertion Site** while allowing airflow for stump desiccation
 - **Easy Repositioning and Removal** of each catheter without removing adhesives
 - **Standardized Method** of securement to save time and simplify the procedure
 - **Sterile Packaging** to ensure a clean environment at the time of securement
- **Simple Application** after catheter insertion to maintain the NICU workflow
 - **Clear Visualization** of catheter position and insertion site for nurses
 - **Soft Medical Grade Silicone** to minimize skin irritation
 - **Full Access** to catheter insertion site for procedures using LifeBubble Buddy
 - **Complete Kit** includes pre-cut neonatal skin-friendly hydrocolloid adhesive

See Instructions for Use for full instructions, warnings, precautions and contraindications.

1. Perl JR, Crabtree-Beach T, Olyaei A, Hedges M, Jordan BK, Scottoline B. Journal of Perinatology 2024.

2. Salimi-Jazi F, Wood L, Rafeeqi T, Yamada N, Fuerch J, Wall J. American Journal of Perinatology 2023.

© 2024 Laborie. All rights reserved.

Ordering information

Description	Product Code	Units Per Box	Contents
LifeBubble Kit - Large	LKL	20	LifeBubble + Adhesive Pack, Large
LifeBubble Kit - Small	LKS	20	LifeBubble + Adhesive Pack, Small
Adhesive Pack - Large	APL	20	Adhesive Pack, Large
Adhesive Pack - Small	APS	20	Adhesive Pack, Small
LifeBubble Buddy	BB	40	LifeBubble Buddy Catheter Insertion Site Access Tool

Laborie¹

USA

T +1 800 522 6743 E usmarketing@laborie.com

F +1 802 878 1122 W laborie.com



Learn more