

EndoTNeedle

Endoscopic Injection Needle

GASTROENTEROLOGY

One hand activation and retraction.
Premium stainless steel construction.
Ultra-sharp needle.



EndoTNeedle

Endoscopic Injection Needle

Introducing EndoTNeedle

Your new everyday needle

EndoTNeedle was engineered to balance performance and value while delivering a highly reliable needle for your most tortuous cases. Designed with quality materials, features from premium needles, and a handle your staff will love. The EndoTNeedle is your new everyday needle.

Click it

- Single-handed use in an ergonomic handle
- Thumb-activated needle extension and spring-loaded retraction

Stick it

- Reliable performance in retroflex and tortuous paths
- Stainless steel sheath prevents kinking
- Smooth round hub which aids in the closure and blockage needed to stop GI bleed cases
- Ultra-sharp needle



Item No.	Minimum Working Channel	Working Length	Needle Extension	Needle Gauge	Units
ET1300	2.8mm	160cm	4mm	25	10
ET1301	2.8mm	230cm	4mm	25	10
ET1302	2.8mm	160cm	5mm	25	10
ET1303	2.8mm	230cm	5mm	25	10
ET1320	2.8mm	160cm	4mm	23	10
ET1321	2.8mm	230cm	4mm	23	10
ET1322	2.8mm	160cm	5mm	23	10
ET1323	2.8mm	230cm	5mm	23	10

Spot is a registered trademark of GI Supply, Inc. EndoT is a trademark of EndoTherapeutics. © 2022 Laborie Medical Technologies Corp. All Rights Reserved.

Laborie¹

USA

T +1 800 522 6743

E gis-gisupplyinfo@laborie.com
W laborie.com



Learn more