

Koala[®] IUPC

Sensor-tipped intrauterine pressure catheter (IUPC)

OBSTETRICS, GYNECOLOGY & NEONATAL
Intrapartum intrauterine monitoring,
optimised for proper placement, easy
true zeroing, small tip, 360° sensing.

Koala Cable Tool

The Koala IUPC is compatible with all major brands of fetal monitors. This tool is designed to help identify what fetal monitor is being used and to ensure the correct cables are ordered for an in-service or evaluation.

The most common cables that are used with the Koala are the **GE Medical (Corometrics) IPC-5014 cable**, the **Philips IPC-5080 cable** and the **Philips IPC-5012 cable**.



▲ IPC-5014



▲ Compatible with **GE Corometrics 250cx**



▲ IPC-5080



FM50



FM30

▲ Compatible with **Philips Avalon Monitors**

As a reminder when zeroing the transducer while using the **Philips Avalon monitor** the catheter needs to be disconnected from the cable, not the cable from the adapter cable.



▲ IPC-5012



▲ Compatible with **All Other Philips Monitors**

Koala Cable Tool



▲ IPC-5016



▲ Compatible with **Spacelabs Medical AMS 94000**



▲ IPC-5090



▲ Compatible with **EDAN Instruments**



The Koala Intrauterine Pressure Catheter is the only IUPC that has no electronics in the catheter.

The external transducer provides easy, true zeroing in utero.

Ordering Information

Product Description	Product Code	Box Quantity	Monitor Cables	Product Code	Box Quantity
Koala Intrauterine Pressure Catheter	IPC-5000E	10/ea	Philips 8040A (40V/V/mmHg)	IPC-5018	1 Cable
			Philips Avalon	IPC-5080	1 Cable
			All other Philips Monitors	IPC-5012	1 Cable
			GE Medical (Corometrics) All Models	IPC-5014	1 Cable
			Air-Shields (Litton) Fetascan 1400	IPC-5016	1 Cable
			Spacelabs Medical (AMS) Spacelabs 94000	IPC-5050	1 Cable
			Neovent Medical Stan S31	IPC-5060	1 Cable
			Edan Instruments, Inc.	IPC-5090	1 Cable

© 2023 Laborie. All rights reserved.

See Instructions for Use for full instructions, warnings, precautions and contraindications

Laborie¹

AUSTRALIA

T +07 3862 2000E orders@laborie.com.au

F +07 3862 2077W aus.laborie.com



Learn more