

RenovaRP[®]

Centesis System

GASTROENTEROLOGY

Safe and effective drainage of abdominal ascites and pleural effusion



RenovaRP[®]
Centesis System


Building on the RenovaRP Legacy

The new RenovaRP Centesis System builds on the legacy of the RenovaRP Paracentesis Pump with the addition of thoracentesis as an intended use.

The new RenovaRP also features a simple, but effective, slide mechanism to retain the RenovaRP Tube Set, which enhances ease of use compared to the previous generation. In addition to the slide, the handle has been redesigned to make the Pump easier to carry and reduce the risk of drops during transport.

The RenovaRP Centesis System is specifically indicated for paracentesis and thoracentesis

Abdominal paracentesis



Thoracentesis



RenovaRP Centesis Pump: Adjustable drainage of abdominal ascites and pleural effusion



- ▲ Engaging the newly designed slider, securing RenovaRP Tubing



- ▲ Newly redesigned RenovaRP Centesis Pump handle

RenovaRP Drainage Bags

The RenovaRP Drainage Bags come individually packaged and sterile, and offer the following benefits when compared to bottles and canisters.

- RenovaRP Drainage Bags have a larger capacity compared to the standard 1 L capacity of many bottles and canisters.
- The small footprint saves valuable storage space and increases portability.

RenovaRP Centesis Kit



The RenovaRP Centesis Kit offers the convenience of a standard procedure kit while including the proprietary RenovaRP Tubing which is required to operate the system.

- The RenovaRP Centesis Kit contains standard items as well as the proprietary RenovaRP Tubing.
- The Tubing features a T-shaped end with dual stopcocks to direct flow from one bag to the other, eliminating the need to start and stop the procedure.
- The Kits are sold in a box of five and come with a 5F Merit One-Step™ Catheter.

RenovaRP Cart

The easily maneuverable RenovaRP Cart with secure pump storage and a place for all your procedural supplies, enhances the portability and convenience for which the RenovaRP Centesis Pump is already known.

- Two full extension drawers and additional pull out work surface.
- Drawers and Pump cabinet lock with one key for added security.
- Ample storage with six side wing shelves, two with dispensing slots.
- Four wheels with two lockable castors.
- 100% total replacement lifetime warranty.



2023 CMS National Payment¹

CPT Code	Description	Physician Office	ASC	Hospital Outpatient
49083	Abdominal Paracentesis with Imaging Guidance	\$300	\$430	\$826 APC 5301
32555	Thoracentesis with Imaging Guidance	\$321	\$301	\$578 APC 5181

RenovaRP Standard Warranty & Service Plans

Coverage	1 Year Standard Warranty	3 Year Service Plan (from time of purchase)	5 Year Service Plan (from time of purchase)
Accessories - Power Cable		✓	✓
Service Loaner		✓	✓
Technical Support		✓	✓
Accidental Damage		✓	✓
Manufacturing Defects	✓	✓	✓
Pricing		\$995	\$1495

Service Plans are available for purchase up to 90 days after purchase of the Renova Pump. Must be used with RenovaRP Tube Set or Service Plan will be void.

To return the unit for service, please contact Laborie Technical Support at service@laborie.com.

RenovaRP Centesis Ordering Information

Item No.	Description	Unit
GIS-35	RenovaRP Centesis Pump	1 Unit
GIS-29	RenovaRP Drainage Bags	Box of 20
GIS-31	RenovaRP Centesis Kit (no catheter)	Box of 5
GIS-32	RenovaRP Centesis Kit (5F catheter)	Box of 5
GIS-34	RenovaRP Cart	1 Unit

1. CPT - Current Procedural Terminology 2022, American Medical Association. Chicago, IL 2022. CPT is a registered trademark of the American Medical Association. One-Step is a trademark of Merit Medical Systems Inc. RenovaRP is a registered trademark of GI Supply, Inc. © 2023 Laborie Medical Technologies Corp. All Rights Reserved.

Laborie¹

USA

T +1 800 522 6743

E gis-orders1@laborie.com

W laborie.com



[Learn more](#)

174096-0 11/23