


# RenovaRP<sup>®</sup> Cart

Part of the RenovaRP Centesis System

GASTROENTEROLOGY  
Convenience  
Portability  
Storage



# Trying to cut down on traffic between departments?

The Renova Cart allows you to securely transport your Renova Pump from department to department with ease so you can treat your patients wherever it is safest to do so.



## RenovaRP cart features

- Two full extension drawers and additional pull out work surface
- Drawers and Pump cabinet lock with one key for added security
- Ample storage with six (6) side wing shelves, two (2) with dispensing slots
- Four (4) wheels with two (2) lockable castors

Upper Drawer  
Width: 11 3/8"  
Depth: 15 3/8"  
Height: 3 1/4"

Lower Drawer  
Width: 11 3/8"  
Depth: 15 3/8"  
Height: 7 3/8"

## Specially designed to simplify Renova Pump storage and transport

The easily maneuverable RenovaRP Cart with secure pump storage and a place for all of your procedural supplies, enhances the portability and convenience for which the RenovaRP Paracentesis Pump is already known.

Item No.	Description	Unit
GIS-34	RenovaRP Cart	1 unit

100% Total Replacement Lifetime Warranty

RenovaRP is a registered trademark of GI Supply, Inc. © 2023 Laborie Medical Technologies Corp. All Rights Reserved.

**Laborie**

USA

T +1 800 500 6743

E gis-orders1@laborie.com  
W laborie.com



Learn more

184189-0